BIP HistoCore®

Biopsy Device, disposable

de	-	Gebrauchsanweisung	5
en	-	Instructions for use	7
fr	-	Mode d'emploi	10
it	-	Istruzioni per l'uso	12
es	-	Instrucciones de uso	15
pt	-	Instruções de utilização	18
nl	-	Gebruiksaanwijzing	20
da	-	Brugsanvisning	23
no	-	Bruksanvisning	25
sv	-	Bruksanvisning	28
fi	-	Käyttöohje	30
et	-	Kasutuskorraldus	33
lv	-	Lietošanas pamācība	35
lt	-	Naudojimo instrukcija	37
pl	-	Instrukcja użytkowania	40
cs	-	Návod na použití	42
sk	-	Návod na použitie	45
sr	-	Uputstvo za upotrebu	47
hu	-	Használati útmutató	50
ro	-	Instrucțiune de utilizare	52
bg	-	Инструкции за употреба	55
el	-	Οδηγίες χρήσης	57
tr	-	Kullanım kılavuzu	60

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ENGLISH

Intended purpose

The biopsy device is intended for the removal of soft tissue samples for histological examination. It is intended for single-use and transient invasive application.

• Indications

Indications are anomalies in soft tissues as breast, prostate, liver, kidney, lung or thyroid.

• Contraindications

It is not to be used in bones, in the cardio vascular system and in the central nervous system.

It is not intended for use on adolescents, children, infants and newborns.

Further contraindications are blood clotting disorder and anticoagulant treatment. In case of increased risks, application depends on decision of healthcare professionals.

• Residual risks and side effects

Residual risks and undesirable side effects are tumor seeding due to displacement of tissues or cells and general puncture-related complications such as bleeding, infection, injury to surrounding tissue/organ and pain as well as pneumothorax and air embolism if used on lung.

• Intended users

Intended users are health care professionals only (physicians).

• Intended patients

Intended patients are adults.

BIP HistoCore®

Biopsy Device, disposable

Clinical benefits

- Better tissue results,
 - by powerful, spring-loaded system, adjustable penetration depth (25 mm standard / 18 mm reduced) and fully- or semi-automatic use for precise positioning.
- Optimized tissue penetration and increased target precision
 - by WingTip-Trocar® tip.
- Better visibility in ultrasound,
 - by echogenic zone at cannula.
- Fewer side effects,
 - by shorter needle tip for reduced trauma to the tissue behind the biopsy area.
- Safe use,
 - by possibility of one-handed operation, indicator of the cocking status and protection against accidental triggering.

Description

The BIP HistoCore® is an automatic, spring-operated disposable biopsy device. It offers a fully-automatic and a semi-automatic mode with two adjustable penetration depths (25 mm / 18 mm) each. The biopsy device is designed for one-handed operation. It consists of a plastic handle with indicator of cocking status and an integrated biopsy needle. The biopsy needle is made of medical-grade stainless steel and consists of an inner stylet with 19 mm sample notch, available with beveled tip (type HC) or with WingTip-Trocar® tip (type HC-T /HC-TX) (see Fig. 3), an outer cannula and an attached protective hose (see Fig. 2). The outer cannula is equipped with centimeter depth marks for direct visual assessment of insertion depth and has an echogenic zone at the distal end for better visibility in ultrasound. The HC-TX type also contains a coaxial cannula pre-mounted on the biopsy needle (see Fig. 2, 10). The BIP HistoCore® is sterile, single-packaged in peel pouches and is available in various diameters and lengths (see Fig. 2, Tab. 1). Respective sales units are additionally packaged in protective cardboard box.

Compatibility

The BIP HistoCore® biopsy device can be used in combination with compatible BIP Coaxial Cannulas (see Tab. 1).
The BIP HistoCore® biopsy device can be used with conventional imaging using ultrasound or X-ray (mammography, computer tomography CT).

Application

1. Preparation of the patient

The procedure is to be performed applying aseptic techniques. The biopsy normally takes place through healthy, uninjured skin. When using a biopsy device with a large diameter, an incision should be made beforehand with a scalpel according to the discretion of the user. In addition to the requisite explanation of risks to the patient, adequate local anesthesia can be made according to the discretion of the user.

2. Remove the biopsy device from the sterile package

Check sterile packaging and the product for expiration date and visual damages before use.
Open the peel pouch at the end with the markings (see Fig. 1):

- First open the seal at the two corners from the centre to the outside 
- Then pull the peel pouch open smoothly  and remove the product aseptically.

3. Cocking of the biopsy device

Hold the device in your hand for this purpose. Using your index finger and middle finger, press the cocking handle back into the housing twice until an audible clicking sound can be heard (Fig. 4).

With the first stroke the cannula slide is cocked and with the second stroke the needle slide is cocked.

NOTE: For application "without pre-shot" the device is cocked in front of the lesion only once.

The current cocking state can be viewed at the indicator window (Fig. 5).

- If no blue bar can be seen, the biopsy device is fully released and not ready for use.
- If one blue bar can be seen in the middle, the cannula slide is cocked, the sample notch is open – the device is ready for use for application "without pre-shot".
- If both blue bars can be seen, the cannula slide and the needle slide are cocked – the device is ready for use in fully- and semi-automatic mode.

4. Setting the penetration depth

The standard setting is 25 mm. In order to select 18 mm, the slide control must be pushed left (Fig. 6).

NOTE: The penetration depth can only be switched while the device is in the fully cocked state. The slide control must always be pushed completely up to the stop.

5. Remove protective hose

Carefully remove protective hose from the needle (see Fig. 2) and place aside.

6. Insertion into tissue

The insertion of the needle into the body should be carried out under the control of a suitable compatible imaging procedure.

When using the biopsy needle without a coaxial cannula, the needle is directly inserted and positioned for tissue removal.

When using the biopsy needle in combination with a coaxial cannula, the coaxial cannula must be inserted first to target and the stylet must be pulled out. The biopsy needle is then inserted through the coaxial cannula until stop and positioned for tissue removal.

When using with integrated coaxial cannula (type HC-TX), the biopsy needle together with the pre-mounted coaxial cannula is inserted to target.

NOTE: The use of BIP coaxial cannula is optional, but recommended in the case of multiple biopsies and in case of application "without pre-shot" in order to minimize trauma and prevent risk of tumor seeding. Make sure that correct coaxial cannula (size and diameter) is used for applied biopsy needle (see Tab. 1, not applicable for type HC-TX).

NOTE: For application "without pre-shot" the stylet is positioned in only once cocked state with the sample notch in the intended target.

7. Performing the biopsy

After the successful positioning of the biopsy needle and without interim change in position, the shot is triggered.

The user can choose between fully-automatic and semi-automatic mode:

BIP HistoCore®

Biopsy Device, disposable

Fully-automatic mode (see Fig. 7):

Press  button, wherein the inner needle and outer cannula are triggered in very short succession and the biopsy process is completely performed.

Semi-automatic mode (see Fig. 7):

Press  button, wherein only the inner needle is shot. The outer cannula is shot with the actuation of the  button. This mode can significantly increase the biopsy yield, depending on the tissue characteristics.

The biopsy needle must never be pulled out of the body without pressing the  button, because a tissue sample could not be removed and there would be the risk of tissue transfer into the punch channel.

Application „without pre-shot“ (see Fig. 7):

For application „without pre-shot“ the device is cocked only once and the cannula is shot after positioning of the needle with the actuation of the  button.

8. Remove biopsy needle

After performing the biopsy procedure, remove the device with the biopsy needle directly or through the coaxial cannula.

NOTE: With BIP HistoCore® type HC-TX, the coaxial cannula must be released after the first biopsy before removing the biopsy needle (Fig. 10). This is achieved by slightly turning the handle of coaxial cannula (Fig. 8) or by pulling back the slider at biopsy device (Fig. 9).

9. Removal of the tissue sample

For the removal of the tissue sample from the sample notch, the outer cannula is retracted by cocking the handle once (Fig. 4). The tissue sample can then be extracted.

If additional biopsies should take place after the removal of the tissue sample, the device is ready for use after re-cocking.

10. Disposal

Dispose of in accordance with local provisions and safety regulations.

Warnings

- For those organs as listed in the indications, safety and performance is demonstrated on sufficient clinical data. If the device is to be used on other similar soft tissues, except those body regions listed as contraindications, safe application depends on decision of healthcare professionals.
- The biopsy device must not be used and must be disposed of if the product (needle and/or housing) is damaged or if the expiration date has elapsed or if the sterile packaging is damaged or unintentionally opened before use.
- Only remove the protective hose immediately prior to performing biopsy to prevent accidental injuries at sharp tip.
- Do not make "empty shots". Activating the biopsy device in the air, not in tissue, can result in damage to the cannula tip and can impair the correct and safe functionality of the biopsy needle.
- •The cocking handle must not be pulled during triggering (shot) to ensure the correct and safe performance.
- Do not apply extreme force to the biopsy needle, especially while stylet is extended from cannula after "single-shot" or for application "without pre-shot". This may cause the stylet or cannula to bend. A strongly bent biopsy needle can impair correct and safe functionality of the product and must be replaced.
- The biopsy device is not MR compatible. For use in MRI environment ensure that the biopsy is only performed outside the 5 Gauss line and MR compatible coaxial cannula (see Tab. 1, CC[xx][yyy]MR) is used for determination of target position.
- The biopsy device is designed for single use and must not be reused. Reuse and re-sterilization impair safe functionality and can lead to cross-infection of patient or material failure of the product.
- The used biopsy device must be disposed of with the protective hose reattached or with other suitable measures to prevent injuries or biohazardous contaminations. The product may only be disposed of fully released.

Precautions

- The warranty applies exclusively to the proper intended use of the product.
- The manufacturer does not recommend any particular biopsy needle size. Healthcare professional users are responsible for selection of proper size depending on patient to be diagnosed.
- Please also note the operating instructions of compatible coaxial cannulas. The manufacturer assumes no liability and warranty for use with non-compatible medical devices.
- Before use inspect each biopsy device for damages at tip or bent cannula that would prevent proper functionality. Do not use damaged or bent needles.
- Use of the biopsy device on dense tissue (e.g. prostate) may cause the needle to bend, especially true for smaller diameter sizes and after performing multiple biopsies. Strongly bent needles must be replaced by new biopsy devices.
- According to the European Medical Device Regulation (2017/745) the manufacturer is obligated to inform you of the following:
All serious incidents that have occurred in relation to the medical device must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.
Note: Serious incidents are understood to mean the death of a patient, user or other person, or the temporary or permanent serious deterioration in the state of health of a patient, user or other person. It does not matter whether these occurred or could occur. The exact definition can be found in MDR 2017/745 Article 2 (65).

Environmental storage conditions

Temperature		10 – 25 °C
Humidity		20 – 50 % r.h. (non-condensing)
Shelf-life		3 years

NOTE:

- Sterility and device functionality are not affected if storage temperature is outside specified environmental conditions over limited short period as long as average storage conditions comply with specifications over respective timeframe.
- For operation no limited environmental conditions apply.
- For transportation packaging can be exposed over limited short period to environmental conditions outside of those as specified for storage.

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©2024 BIP GmbH, Germany, all rights reserved.