



Silicone Carving Blocks

Implantech offers implantable silicone carving blocks which allow the surgeon to easily carve customized implants for their patient's specific needs. Our blocks are available in four different styles, including rectangular and preformed, and may be ordered in various sizes, durometers and colors.

Uses include aesthetic enhancement, reconstruction following trauma or correcting deficient areas of tissue due to congenital defects.

Implantech[®]


Superior Patient Aesthetics

800.733.0833 | implantech.com

Silicone Carving Blocks


Contoured Carving Block - Style 1

CATALOG NUMBER	DIMENSIONS (in centimeters)			DUROMETER	COLOR	VOLUME
A	B	C				
CCB1-1	11.9	5.3	1.5	Extra Soft	Clear	43cc
CCB1-2	12.6	6.1	1.6	Extra Soft	Clear	58cc
CCB1-3	13.7	6.8	1.7	Extra Soft	Clear	73cc
CCB1-4	14.0	7.8	1.9	Extra Soft	Clear	102cc
CCB1-5	16.2	9.4	2.0	Extra Soft	Clear	149cc




Contoured Carving Block - Style 3

CATALOG NUMBER	DIMENSIONS (in centimeters)			DUROMETER	COLOR	VOLUME
A	B	C				
CCB3-1	11.5	4.7	2.1	Extra Soft	Clear	55cc
CCB3-2	14.1	5.0	1.8	Extra Soft	Clear	58cc
CCB3-3	17.2	6.1	2.1	Extra Soft	Clear	93cc
CCB3-4	20.3	6.5	2.1	Extra Soft	Clear	130cc




Contoured Carving Block - Style 4

CATALOG NUMBER	DIMENSIONS (in centimeters)			DUROMETER	COLOR	VOLUME
A	B	C				
CCB4-1	2.60	2.0	-	Extra Soft	Clear	5cc
CCB4-2	3.25	2.5	-	Extra Soft	Clear	10cc
CCB4-3	3.90	3.0	-	Extra Soft	Clear	19cc
CCB4-4	4.55	3.5	-	Extra Soft	Clear	29cc
CCB4-5	5.00	4.0	-	Extra Soft	Clear	40cc




Contoured Carving Block - Style 8

CATALOG NUMBER	DIMENSIONS (in centimeters)			DUROMETER	COLOR	VOLUME
A	B	C				
CCB8-1	2.60	2.0	-	Extra Extra Soft	Clear	5cc
CCB8-2	3.25	2.5	-	Extra Extra Soft	Clear	10cc
CCB8-3	3.90	3.0	-	Extra Extra Soft	Clear	19cc
CCB8-4	4.55	3.5	-	Extra Extra Soft	Clear	29cc
CCB8-5	5.00	4.0	-	Extra Extra Soft	Clear	40cc



Rectangular Carving Block

CATALOG NUMBER	DIMENSIONS (in centimeters)			DUROMETER	COLOR	VOLUME
A	B	C				
BL-1.0-NS	13.6	5.8	1.0	*	*	-
BL-1.3-NS	14.8	7.6	1.3	*	*	-
BL-2.0-NS	14.9	7.7	2.0	*	*	-
BLK-10-NS	12.0	7.0	1.5	Soft	Clear	-
BLK-20-NS	12.0	7.0	1.5	Medium	Clear	-
BLK-40-NS	12.0	7.0	1.5	Firm	Clear	-



*Specify hardness and color when ordering rectangular blocks.
Available hardness: Very Soft, Soft, Medium, Firm, and Very Firm.
Available colors: Clear, Opaque, and Fleshtone.
May take 2-4 weeks processing time.

NOTE: All dimensions and volumes are approximate.
Implants may be trimmed to meet the individual patient needs.

Implantech®
Superior Patient Aesthetics

6025 Nicolle Street • Ventura, CA 93003

800.733.0833 | implantech.com

Individual sizes for all contoured carving blocks available (Non-Sterile)

ABC 007/1120