

UROMED transurethraler Silikonballonkatheter »soft«

UROMED Transurethral Silicone Balloon Catheter »soft«

REF 1095

Nelatonspitze
2 Augen


Nelaton tip
2 eyes

Silikonballonkatheter zur transurethralen Blasendrainage

- kurze, konische Spitze
- großlumig
- Soft-Ballon
- reines Silikon
- mit Katheterstopfen
- sehr weich

Silicone balloon catheter for transurethral bladder drainage

- short, conical tip
- wide lumen
- soft balloon
- pure silicone
- with catheter plug
- very soft



REF REF	Ch./Fr. Ch./Fr.	Ø Ø	Ballon Balloon	Länge Length	Effektivlänge Effective length	PZN CPN
1095.12	12	4,0 mm	3-6 ml	430 mm	369 mm	04741526
1095.14	14	4,7 mm	3-8 ml	430 mm	369 mm	04741532
1095.16	16	5,3 mm	5-10 ml	430 mm	369 mm	04741549
1095.18	18	6,0 mm	5-10 ml	430 mm	369 mm	04741555
1095.20	20	6,7 mm	5-10 ml	430 mm	369 mm	04741561
1095.22	22	7,3 mm	5-10 ml	430 mm	369 mm	04741578

Hilfsmittelkatalog: 15.25.15.6030

Verpackung: steril
Originalkarton: 10 Stück

Packaging: sterile
Original carton: 10 pieces