



ITABLOC



Boccaglio per gastroenterologia

BITEBLOCK FOR GASTROENTEROLOGY

Meditalia®



Meditalia®

BOCCAGLI PER GASTROENTEROLOGIA

BITEBLOCKS FOR GASTROENTEROLOGY

Cintino premontato ad elasticità differenziata
Pre-assembled headstrap with differentiated elasticity

Struttura rigida
Rigid construction

Abbassalingua munito di grip
Tongue pad with grip

Ampi accessi laterali
Large side access

Stampo cromato garantisce angolature smusse
Chromed mold guarantees smooth angles

Products info

References	Oral diameter	Material	Colour	Headstrap
MED/102	21x25 mm	polypropylene	yellow	•
MED/102/B	21x25 mm	polypropylene	yellow	
MED/102/A	21x25 mm	polycarbonate 	transparent	•
MED/102/A/B	21x25 mm	polycarbonate 	transparent	
MED/106	18x20 mm	polypropylene	purple	•
MED/106/B	18x20 mm	polypropylene	purple	
MED/106/P	16x20 mm	polypropylene	orange	•
MED/106/P/A	16x20 mm	polycarbonate 	transparent	•