

BASKET DA RECUPERO 4/8 FILI

RETRIEVAL BASKET 4/8 WIRES

DESTINAZIONE D'USO

Cestelli utilizzati per la rimozione di polipi e corpi estranei durante procedure di endoscopia.

CARATTERISTICHE

Versione rotante 360°. Punta a pallina atraumatica. Ottima forza radiale. Sicurezza e memoria di forma. Catetere in acciaio inox intrecciato che garantisce un'ottima flessibilità. Morfologia a 4 fili. Testa a rete a 8 fili che permette una sicura e veloce presa. Versione a 8 fili. Confezione sterile monouso. Box da 10 pz.

INDICATIONS

Basket Used for retrieval of polyps and foreign bodies during endoscopy procedures.

FEATURES

360° Rotatable version. Atraumatic ball-tip. Excellent radial force. Security and memory shape. Catheter in braided stainless steel that guarantees excellent flexibility. 4 wire basket with 8 wire head-net for a safe and quick retrieval. 8 wire version. Packaged sterile for, single use. 10 pcs per box.

PRODUCT INFO

References	Working length	Sheath Diameter	Basket	Diameter	Channel
MED-145-BAS	215 cm	2,3 mm	4/8 wire hexag	20 mm	2,8 mm
MED-146-BAS	215 cm	2,3 mm	4/8 wire hexag	30 mm	2,8 mm