


PICA

Pigtail Universal Drainage Catheter

- Used for nephrostomy, abscess, and other multipurpose drainage
- Catheter placement is performed using one step (direct access) technique.
- Hydrophilic coating to reduce friction during insertion
- Available with StringLock
- Large oval holes in pigtail for maximum drainage




HYDROPHILIC DRAINAGE CATHETER SET-DIRECT ACCESS Without Safety String Lock (Box of 5)

Catheter French Ø	Order Code	Accepts Guide Wire Ø	Catheter Length
6F	100106YY	0.035" (0.89mm)	25cm (YY=25) 30cm (YY=30)
7F	10010730	0.035" (0.89mm)	
8F	100108XX	0.038" (0.97mm)	
10F	100110XX	0.038" (0.97mm)	
12F	100112YY	0.038" (0.97mm)	30cm
14F	10011430	0.038" (0.97mm)	

To Order Single Pigtail Drainage Catheter use following code:
1001FFLL (FF=French and LL=Length)

SET INCLUDES:

- Pigtail Drainage Catheter
- Metal Stiffening Cannula
- Trocar Stylet
- Split Straightener
- ON REQUEST AVAILABLE ALSO WITH CONNECTING TUBE



HYDROPHILIC DRAINAGE CATHETER SET-DIRECT ACCESS With Safety String Lock (Box of 5)

Catheter French Ø	Order Code	Accepts Guide Wire Ø	Catheter Length
6F	151106YY	0.035" (0.89mm)	25cm (YY=25) 30cm (YY=30)
7F	15110730	0.035" (0.89mm)	
8F	151108YY	0.038" (0.97mm)	
10F	151110YY	0.038" (0.97mm)	
12F	151112YY	0.038" (0.97mm)	30cm
14F	15111430	0.038" (0.97mm)	

To Order Single Pigtail Drainage Catheter use following code:
1511FFLL (FF=French and LL=Length)

SET INCLUDES:

- Sting Lock Pigtail Drainage Catheter
- Metal Stiffening Cannula
- Trocar Stylet
- Split Straightener
- ON REQUEST AVAILABLE ALSO WITH CONNECTING TUBE

