



Hybrida[®]



Ansa da polipectomia

POLYPECTOMY SNARE

Meditalia[®]



Meditalia

ANSA DA POLIPECTOMIA PLURIFILAMENTO

BRAIDED POLIPECTOMY SNARE

Tipologia con rotazione graduata 1:1, tramite filo in nitinol inserito all'interno del catetere

Graduated rotation 1:1 version, through a nitinol wire inside the catheter

I fili intrecciati sono di spessore maggiorato rispetto agli standard di mercato, ciò determina una migliore performance nel taglio e memoria di forma del loop

The greater thickness of the braided wires compared to the market standards allows a better cutting performance and a better memory shape of the loop

Gamma completa di misure da 7 a 50 mm nelle morfologie ovale, esagonale e crescent

Full range of sizes from 7 to 50 mm with oval, hexagonal and crescent shape

Products info

References	Working length	Sheath diameter	Channel	Loop diameter	Loop shape
MED-409-SNA	230 cm	2,3 mm	2,8 mm	7 mm	oval
MED-180-SNA	230 cm	2,3 mm	2,8 mm	10 mm	oval
MED-120-SNA	230 cm	2,3 mm	2,8 mm	20 mm	oval
MED-121-SNA	230 cm	2,3 mm	2,8 mm	30 mm	oval
MED-122-SNA	230 cm	2,3 mm	2,8 mm	50 mm	oval
MED-181-SNA	230 cm	2,3 mm	2,8 mm	10 mm	hexagonal
MED-123-SNA	230 cm	2,3 mm	2,8 mm	20 mm	hexagonal
MED-124-SNA	230 cm	2,3 mm	2,8 mm	30 mm	hexagonal
MED-182-SNA	230 cm	2,3 mm	2,8 mm	10 mm	crescent
MED-126-SNA	230 cm	2,3 mm	2,8 mm	20 mm	crescent
MED-127-SNA	230 cm	2,3 mm	2,8 mm	30 mm	crescent
MED-208-SNA	230 cm	2,3 mm	2,8 mm	20 mm	oval 
MED-209-SNA	230 cm	2,3 mm	2,8 mm	30 mm	oval 
MED-210-SNA	230 cm	2,3 mm	2,8 mm	20 mm	hexagonal 
MED-211-SNA	230 cm	2,3 mm	2,8 mm	20 mm	crescent 