

UROMED transurethraler Silikonballonkatheter für Frauen

UROMED Transurethral Silicone Balloon Catheter for Women

REF 1370

Nelatonspitze
2 Augen

Nelaton tip
2 eyes

Silikonballonkatheter zur transurethralen Blasendrainage

- kurze, konische Spitze
- großlumig
- Spitze röntgenologisch darstellbar
- Soft-Ballon
- reines Silikon

Silicone balloon catheter for transurethral bladder drainage

- short, conical tip
- wide lumen
- radiopaque tip
- soft balloon
- pure silicone



REF REF	Ch./Fr. Ch./Fr.	Ø Ø	Ballon Balloon	Länge Length	Effektivlänge Effective length	PZN CPN
1370.12	12	4,0 mm	3-8 ml	223 mm	154 mm	03280862
1370.14	14	4,7 mm	3-8 ml	224 mm	157 mm	03280879
1370.16	16	5,3 mm	5-10 ml	225 mm	158 mm	03280885
1370.18	18	6,0 mm	5-10 ml	225 mm	158 mm	03280891
1370.20	20	6,7 mm	5-10 ml	226 mm	159 mm	03280916
1370.22	22	7,3 mm	5-10 ml	226 mm	159 mm	03280922
1370.24	24	8,0 mm	5-10 ml	227 mm	160 mm	03280939

Hilfsmittelkatalog: 15.25.15.6086

Verpackung: steril
Originalkarton: 10 Stück

Packaging: sterile
Original carton: 10 pieces