

EndoView CY Series Single use Flexible Cystoscope

Technical Specifications

Model		CY-M32	CY-M40	CY-M50	CY-M52
Optical system	Field view	120°			
	Depth of field	3-50 mm			
Insertion tube	Distal end outer diameter	3.2mm	4.0mm	5.0mm	5.2mm
	Working length	350 mm			
Working Channel		1.2mm	2.0mm	2.8mm	2.2mm
Angulation	Bending	210 ° U/D			
Total Length		605mm			

Processor	
Function	Photo
	Video
	White Balance
	Freeze
Storage	Minimal 16G
Connection	DVI, HDMI
	VGA
	SDI, AV
Display	
Screen	LCD
	Size: 8 inch
Function	Photo
	Video
	White Balance
	Freeze
Battery	Type: Lithium
	Maximum 4 hours
Connection	HDMI, DVI

Single use Flexible Cystoscope	
Order Information	Product
2017B-PA004 EU	CY-M32
2017B-PA005 US	
2017B-PA003 EU	CY-M40
2017B-PA006 US	
2017B-PA001 EU	CY-M50
2017B-PA007 US	
2017B-PA002 EU	CY-M52

Shenzhen HugeMed Medical Technical Development Co.,Ltd

Add: 416-1,516-1, Building 2, No. 1, Mawu Road, Baoan
Community, Yuanshan Street, Longgang District,
Shenzhen 518115, China
Tel: +86-755-2227 5833
Fax: +86-755-2227 5833
Email: sale@hugemed.net

www.hugemed.net