

# UROMED NEPHROSTOMIE BALLON-KATHETER

## REF 2525 / 2526

Zentral offene Spitze  
2 Augen, Länge 43 cm

### Ballon-Katheter zur perkutanen Nierendrainage


- kurze, atraumatische Spitze mit abgerundeter Kante
- flache Ballonform
- röntgenologisch darstellbar
- graduert
- Soft-Ballon
- reines Silikon
- halbfest

### Balloon Catheter for percutaneous kidney drainage

- short, atraumatic tip with rounded edge, central open, 2 eyes, length 43 cm
- flat shape of balloon
- overall radiopaque
- graduated
- soft balloon
- pure silicone
- semi-solid

Zur sicheren Blockung ohne Diffusion  
verwenden Sie bitte unser Füllmedium  
REF 1550.05.

*For a safe inflation without  
diffusion please use our  
filling medium REF 1550.05.*



REF 2525

REF 2526

2 zusätzliche Löcher zur Kelchdrainage  
2 additional holes for calyceal drainage

| REF<br>REF | Größen<br>Sizes | Ballon<br>Balloon | PZN<br>PZN |
|------------|-----------------|-------------------|------------|
| 2525.10    | Ch. 10          | 2 ml              | 03281057   |
| 2525.12    | Ch. 12          | 5 ml              | 03281063   |
| 2525.14    | Ch. 14          | 5 ml              | 03281086   |
| 2525.16    | Ch. 16          | 5 ml              | 03281092   |
| 2525.18    | Ch. 18          | 5 ml              | 03281100   |
| 2525.20    | Ch. 20          | 5 ml              | 03281117   |
| 2525.22    | Ch. 22          | 5 ml              | 03281123   |
| 2525.24    | Ch. 24          | 5 ml              | 03281146   |

| REF<br>REF | Größen<br>Sizes | Ballon<br>Balloon | PZN<br>PZN |
|------------|-----------------|-------------------|------------|
| 2526.10    | Ch. 10          | 2 ml              | 04741868   |
| 2526.12    | Ch. 12          | 5 ml              | 03281152   |
| 2526.14    | Ch. 14          | 5 ml              | 03281169   |
| 2526.16    | Ch. 16          | 5 ml              | 03281175   |
| 2526.18    | Ch. 18          | 5 ml              | 03281181   |
| 2526.20    | Ch. 20          | 5 ml              | 03281198   |
| 2526.22    | Ch. 22          | 5 ml              | 03281206   |
| 2526.24    | Ch. 24          | 5 ml              | 03281212   |

Verpackung: doppelt, steril  
Original-Karton: 10 Stück

Packing: double, sterile  
Original carton: 10 pieces