

UROMED BALLON-KATHETER-SET »NEPHROquick®«

REF 2608/ 2530

Pigtail war gestern!


Pigtail – Solution of the Past!

Ballon-Katheter-Set zur perkutanen Nierendrainage

- Dilator und Ballon-Katheter in einem
- Ballon-Katheter aus reinem Silikon
- röntgenologisch darstellbar
- graduert
- knickstabil
- halbfest

Balloon Catheter Set for percutaneous kidney drainage

- dilator and balloon catheter in one
- balloon catheter made of pure silicone
- overall radiopaque
- graduated
- nonkinking
- semi-solid



REF 2700

Lunderquist
Führungsdraht
0,035"/0,89 mm
Länge 90 cm mit
Einführhilfe
Lunderquist
guidewire
0,035"/0,89 mm
length 90 cm with
insertion guide

REF 2613

Punktionskanüle
17,5 GA / 1,3 mm
Stichlänge
18,5 cm
puncture
cannula
17,5 GA / 1,3 mm
needle length
18,5 cm

REF 2650

Fascial-Dilatator
Ch. 8 / 2,7 mm
fascial dilator
Ch. 8 / 2,7 mm

REF 2650

Fascial-Dilatator
Ch. 10 / 3,3 mm
fascial dilator
Ch. 10 / 3,3 mm

REF 2530

»NEPHROquick®« Katheter und Obturator mit
Katheter-Konnektor, »INTEGRAL« Ch. 10,
Ballon 2 ml
»NEPHROquick®« Catheter and obturator with
catheter connection, »INTEGRAL« Ch. 10,
balloon 2 ml

Alle Nephrostomie-Ballon-Katheter-Sets
enthalten außerdem das UROMED-
Füllmedium zum Blocken
des Ballons.

In addition the UROMED
Filling Medium for inflating the balloon is
included in all nephrostomy balloon catheter sets.



Das erste Nephrostomie-Set mit
einem 10 Charrière Ballon-Katheter.
Es lässt sich in nur wenigen Schritten
applizieren, vereint aber dennoch
alle Vorteile der bisher bekannten
Techniken:

- sichere Fixierung und geringste
Raumforderung im Nierenbecken
 - hoher Tragekomfort
 - gute Dilatation
- Überzeugen Sie sich von
»NEPHROquick®«.

The first nephrostomy set with a
10 Charrière balloon catheter.
It can be applied with a few steps
only, however, it combines all
advantages of the techniques
known so far:

- safe fixation and lowest space
requirement in the renal pelvis
 - high wearing comfort
 - good dilation
- Go and find out yourself about
»NEPHROquick®«.

REF REF	Größe Size
2608.10	Ch. 10
2530.10	Ch. 10

Verpackung: doppelt, steril
Original-Karton: 5 Stück (REF 2608)
10 Stück (REF 2530)

Packing: double, sterile
Original carton: 5 pieces (REF 2608)
10 pieces (REF 2530)